#### WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



### INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification 6:

(11) International Publication Number:

WO 98/37949

B01D 61/00, 63/00

(43) International Publication Date:

3 September 1998 (03.09.98)

(21) International Application Number:

PCT/US98/03696

A1

(22) International Filing Date:

25 February 1998 (25.02.98)

(30) Priority Data:

60/038,909

26 February 1997 (26.02.97)

US

(71) Applicant (for all designated States except US): MILLIPORE CORPORATION [US/US]; 80 Ashby Road, Bedford, MA 01730 (US).

(72) Inventors; and

- (75) Inventors/Applicants (for US only): KOPACIEWICZ, William [US/US]; 32 Archalaus Place, West Newbury, MA 01985 (US). SHEER, Donald, G. [US/US]; 3 Hemlock Lane, Norfolk, MA 02056 (US). ARNOLD, Todd, E. [US/US]; 28 Columbus Avenue, Beverly, MA 01915 (US). GOEL, Vinay [US/US]; 9 Washington Drive, Acton, MA 01720 (US).
- (74) Agents: LEMACK, Kevin, S. et al.; Nields, Lemack & Dingman, Suite 8, 176 E. Main Street, Westboro, MA 01581 (US).

(81) Designated States: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, GW, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO,-RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG).


**Published** 

With international search report.

(54) Title: CAST MEMBRANE STRUCTURES FOR SAMPLE PREPARATION

#### (57) Abstract

A method for casting-in-place composite and/or nonfilled structures (11) which are useful as sorptive or reactive media or for size-based separations. Any particular housing size or configurationcan be used, and the inclusion of a large amount of adsorptive particles in polymer is achieved while still maintaining the membrane three dimensional structure. In a first preferred embodiment, the composite structures comprise particles entrapped within a porous polymeric substrate, and are cast-in-place into a housing such as a pipette tip, thereby providing an effective platform for micromass handling. With the apropriate selection of particle chemistry, virtually any separation or purification operation can be conducted, including selective bind/elute chromatography operations, on sample mass loads less than 1 microgram in volumes. The present invention also encompasses the composite structures as well as sample preparation devices containing the same. In a second preferred embodiment, self-retaining, self-supporting structures are cast in situ in a suitable housing (12) and can be used for size-based separations wherein the cast structure (11) acts as a semipermeable membrane barrier. The present invention also encompasses these structures as well as housings containing the structures.



#### FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

4.7	A Marain	ES	Spain	LS	Lesotho	SI	Siovenia
AL	Albania	FI	Finland	LT	Lithuania	SK	Slovakia
AM	Armenia	FR	France	LU	Luxembourg	SN	Senegal
AТ	Austria		Gabon	LV	Latvia	SZ	Swaziland
AU	Australia	GA		MC	Monaco	TD	Chad
AZ	Azerbaijan	GB	United Kingdom	MD	Republic of Moldova	TG	Togo
BA	Bosnia and Herzegovina	GK	Georgia	MG	Madagascar	TJ	Tajikistan
BB	Barbados	GH	Ghana		The former Yugoslav	TM	Turkmenistan
BE	Belgium	GN	Guinea	MK		TR	Turkey
BF	Burkina Faso	GR	Greece		Republic of Macedonia	TT	Trinidad and Tobago
BG	Bulgaria	HU	Hungary	ML	Mali	UA .	Ukraine
BJ	Benin	1E	<b>Freland</b>	MN	Mongolia		<del>-</del>
BR	Brazil	ΙL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KB	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	ZW	Zimbabwc
CI	Côte d'Ivoire	KР	Democratic People's	NZ	New Zealand		
_	Cameroon		Republic of Korea	PL	Poland		
CM	China	KR	Republic of Korea	PT	Portugal		
CN		KZ	Kazakstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic		Liechtenstein	SD	Sudan		
DE	Germany	LI		SE	Sweden		
DK	Denmark	LK	Sri Lanka	SG	Singapore		
EE	Estonia	LR	Liberia	50	Singapore		
							•

# CAST MEMBRANE STRUCTURES FOR SAMPLE PREPARATION BACKGROUND OF THE INVENTION

5

10

15

20

25

A number of analytical procedures have been developed in the biochemical art wherein it is required to remove solvent from peptide solutions in order to have a more concentrated peptide sample which can be analyzed effectively, or to remove low molecular weight ions or solutes. Many other analytical procedures, involving not only peptides but macromolecular species in general, also have been developed wherein it is necessary to concentrate and/or "desalt" a macromolecular component in a liquid sample, as there is commonly a need in biochemistry/medicinal chemistry for pure analytes devoid of salts, detergents and other contaminants. The presence of these substances can be deleterious, in that they often interfere with subsequent chemical analyses. Analogous situations exist in the environmental art and in chemical analysis.

U.S. Patent No. 4,755,301 discloses a centrifugal method and apparatus for concentrating macromolecules without filtering to dryness. A semipermeable ultrafiltration membrane separates a sample reservoir from a filtrate cup, and filtrate ducts below the membrane are offset sufficiently inward from the edge of the membrane so that when the apparatus is used in a fixed angle centrifuge rotor, filtration stops once the retentate meniscus reaches the centrifugal radial level of the outermost edge of the outermost filtrate duct.

Such ultrafiltration devices are commonly used for the "purification" and/or sample preparation of biomolecules and natural products. For such a process to be successful, a

membrane must be selected that retains the molecules of interest, yet passes the impurities. Although this scenario is relatively straightforward for analytes greater than about 10,000 molecular weight, it becomes increasingly problematic for substances less than about 5000 molecular weight. The reason is due to the fact that the required membrane porosity to retain the about 5000 molecular weight analyte is so low that the water permeability (flow rate) becomes poor and processing times too long. For example, a typical centrifugal "spin time" for a device using a membrane suitable for analytes having a molecular weight of 30,000 or more is about one hour, whereas as many as six hours may be required for analytes of about 1000 molecular weight. Furthermore, such long term exposure to high g-forces frequently results in device failure.

5

10

15

20

25

The sample quantities now common in the art are in the 0.01 to 10 microgram range. At such low loads, efficient sample handling is crucial to avoid loss. Conventional methods and devices for sample preparation are not practical for handling the "microseparation" of such small sample volumes. In addition, ultrafiltration can only effectively concentrate and desalt, and thus the application of adsorption technology at this scale could offer an entirely new approach to micro-mass sample preparation.

One conventional method for making sample preparation devices is to first insert a precut porous plug obtained from, for example, a fiberous glass or cellulose sheet, into the tip of a pipette, followed by the addition of loose particles and a second porous plug, as illustrated in Figure 1. The plugs

serve to retain the particles in place in the pipette tip. However, the plugs also entrap excess liquid thereby creating dead space or volume (i.e., space not occupied by media or polymer that can lead to poor sample recovery, contamination such as by sample carry-over, etc.). However, these procedures cannot be used with extremely small liquid delivery devices such as pipette tips, as there is no practical way to load either the plug or the particles to obtain a microadsorptive device that contains 10 milligrams or less of adsorbent to be used for the aforementioned extremely small sample loads.

5

10

15

20

25

Alternatively, a micro sample preparation device can be made by lodging media in a capillary pipette. However, the flow through such devices is typically slow.

Moreover, although from a mass adsorption standpoint, adsorptive powders offer the highest capacity, they are difficult or indeed impossible to handle in milligram quantities. Although polymer-based adsorptive membrane sheets are relatively easy to handle, their capacity is poor as a result of relatively low substructure surface area.

It is therefore an object of the present invention to provide a sample preparation device which can concentrate, purify and/or desalt molecules from sample solutions.

It is another object of the present invention to provide a sample preparation device which can concentrate, purify and/or desalt molecules from very small sample solutions.

It is another object of the present invention to provide a sample preparation device which can concentrate, purify and/or desalt molecules from sample solutions in a variety of form

geometries.

5

10

15

20

25

It is a further object of the present invention to provide a sample preparation device which can concentrate, purify and/or desalt molecules from very small sample solutions in a variety of form geometries.

It is a still further object of the present invention to provide a sample preparation device that is simple and economic to manufacture.

It is yet a further object of the present invention to provide a method of casting particles in a housing in a variety of housing sizes or geometries.

It is a further object of the present invention to provide a castable membrane that assumes the shape of the housing in which it is cast, and can be retained in that housing without the use of porous plugs.

It is another object of the present invention to provide a castable membrane on a support or substrate.

#### SUMMARY OF THE INVENTION

The problems of the prior art have been overcome by the present invention, which provides a method for casting in-place composite (filled) and/or non-filled structures which are useful as sorptive or reactive media or for size-based separations. In one embodiment, the structures are monolithic and/or continuous. The invention is applicable to a variety of particular housing sizes and configurations, and provides a means of affixing media in a housing of a variety of volumes. The invention enables the inclusion of a substantial (relative

to the increase in surface area of the precipitated structure) amount of media in the polymer while still retaining a three dimensional polymeric structure.

5

10

15

20

25

In a first preferred embodiment, the composite structures particles entrapped within porous polymeric a substrate, such as that shown in Figure 2B, and are cast inplace into a housing of a variety of sizes, such as a pipette tip as illustrated in Figure 2A, thereby providing an effective platform for micromass handling. With the appropriate selection of particle chemistry, virtually any separation or purification operation can be conducted, including selective bind/elute chromatography operations, on sample mass loads less than 1 microgram in volumes of a few microliters, as well as larger The present invention also encompasses mass loads and volumes. the composite structures as well as sample preparation devices containing the same.

In a second preferred embodiment, unfilled structures which may be self-retaining and/or self-supporting are cast in situ in a suitable housing and can be used for size-based separations wherein the cast structure acts as a semi-permeable barrier, or for adsorption. The present invention also encompasses these structures as well as housings containing these structures, such as that shown in Figure 3. The device in Figure 3 is a centrifugal device having a reservoir, a base, and a porous fabric sealed between the reservoir and base. The structures of the invention are cast-in-place on the porous fabric. The device is placed in a suitable vial during operation, and the flux of the device is driven by centrifugal force.

### BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 is a schematic diagram of an adsorptive pipette tip assembled with particles between two porous plugs;

Figure 2A is a schematic diagram of an adsorptive pipette tip prepared by a cast-in-place process in accordance with the present invention;

5

10

15

20

25

Figure 2B is a scanning electron micrograph of a particle loaded cast-in-place structure;

Figure 3 is an further embodiment of a housing into which a cast-in-place structure is added;

Figure 4 is a scanning electron micrograph of a particle loaded membrane prepared with spherical silica gel and polysulfone binder;

Figure 5A illustrates a multi-well array as a suitable housing for the cast-in-place structures of the present invention;

Figure 5B is a side view of an underdrain that can be used with the multi-well array of Figure 5A;

Figure 5C is a top view of a single well of the multi-well array of Figure 5A;

Figure 5D is a cross-sectional view of a well of the underdrain of Figure 5B;

Figure 6 is a reversed phase chromatogram of a 5 peptide mixture;

Figure 7 is a reversed phase chromatogram of a 5 peptide mixture which was bound and eluted from a pipette tip containing cast-in-place C18 silica;

Figure 8 is a reversed phase chromatogram of a 5 peptide mixture which was bound and eluted from a pipette tip containing cast-in-place styrene sulfonate coated silica;

Figure 9 is a reversed phase chromatogram of cytochrome c;
Figure 10 is a reversed phase chromatogram of cytochrome
c after a 15 minute exposure to a pipette tip containing castin-place immobilized trypsin beads;

5

10

15

20

25

Figure 11 is an electrophoretic gel of a rabbit immunoglobulin-g sample that has been bound and eluted from a pipette tip containing cast-in-place immobilized Protein A beads;

Figure 12 is an electrophoretic gel of supercoiled plasmid DNA which has been bound and eluted from a 1000  $\mu$ l pipette tip containing cast-in-place silica;

Figure 13 is an electrophoretic gel of linear DNA fragments ranging from 100-2000 bp which have been bound and eluted from a 200  $\mu$ l pipette tip containing cast-in-place silica;

Figure 14 is an electrophoretic gel of PCR amplified DNA which has been bound and eluted from a 200  $\mu l$  pipette tip containing cast-in-place fumed silica;

Figure 15 is an electrophoretic gel of linear DNA fragments which have been bound and eluted from a 200  $\mu$ l pipette tip containing loose silica and a cast-in-place membrane barrier;

Figure 16A is a scanning electron micrograph of a longitudinal section of a pipette tip containing a cast-in-place particle loaded membrane prepared with spherical silica gel and polysulfone binder; Figure 16B is a scanning electron micrograph of a cross-section of a pipette tip containing a

cast-in-place particle loaded membrane prepared with spherical silica gel and polysulfone binder; and

Figure 17 is a longitudinal section of a pipette tip containing a cast-in-place unfilled uncut membrane.

5

10

15

20

25

#### DETAILED DESCRIPTION OF THE INVENTION

The term "membrane" as used herein includes permeable and semi-permeable three dimensional structures with or without particles, having a porosity suitable for the desired application. The term "composite structure" as used herein includes filled membranes.

In the first preferred embodiment of the present invention, those skilled in the art will recognize that many different particles can be used in the composite structures, depending upon the desired objectives of the resulting device. In the case of adsorptive devices, the ideal device will have rapid adsorption kinetics, a capacity and selectivity commensurate with the application, and allows for elution of bound analyte with an appropriate desorption agent. Suitable adsorptive are polymer bound, particle structures composite adsorptive membrane structures, such as those comprised of chromatographic beads which have been adhered together with a A suitable polymer bound particle laden adsorptive membrane is illustrated in Figure 4. This membrane is comprised of about 80% w/w silica and 20% w/w polysulfone binder, and is produced by Millipore Corporation. A similar membrane is shown in Figure 16A cast-in-place in a pipette tip 50.

5

10

15

20

25

composite structures comprising other micron-size (e.g., 1-30 microns) resin particles derivatized with other functional groups are also beneficial, including styrenedivinyl-benzenebased media (unodified or derivatized with e.g., sulphonic acids, quaternary amines, etc.); silica-based media (unmodified or derivatized with  $C_2$ ,  $C_4$ ,  $C_6$ ,  $C_8$ , or  $C_{18}$  or ion exchange functionalities), to accommodate a variety of applications for peptides, proteins, nucleic acids, and other organic compounds. Those skilled in the art will recognize that other matrices with alternative selectivities (e.g., hydrophobic interaction, affinity, etc.) can also be used, especially for classes of The term "particles" as used molecules other than peptides. herein is intended to encompass particles having regular (e.g., spherical) or irregular shapes, as well as shards, fibers and including metal powders, plastic powders powdered polystyrene), normal phase silica, fumed silica and For example, the addition of fumed silica activated carbon. into a polysulfone polymer results in increased active surface area and is suitable for various applications. Polysulfone sold under the name UDEL P3500 and P1700 by Amoco is particularly preferred in view of the extent of the adherence of the resulting composite structure to polyolefin housing, including polypropylene, polyethylene and mixtures thereof. Other suitable polymer binders include polyethersulfone, cellulose acetate, cellulose acetate butyrate, acrylonitrile PVC copolymer (sold commercially under the name "DYNEL"), polyvinylidene fluoride (PVDF, sold commercially under the name "KYNAR"), olystyrene and polystyrene/acrylonitrile copolymer,

Adhesion to the housing can be enhanced or an analogous effect achieved with these composite structures by means known to those skilled in the art, including etching of the housing, such as with plasma treatment or chemical oxidation; mechanical aids such as rims inside the housing; and inclusion of additives into the housing material that promote such adhesion. Adhesion allows uniform precipitation during casting.

5

10

15

20

25

Devices in accordance with the present invention may incorporate a plurality of composite structures having resin materials with different functional groups to fractionate analytes that vary by charge, size, affinity and/or hydrophobicity; alternately, a plurality of devices containing different individual functional membranes may be used in combination to achieve a similar result. Similarly, one or more membranes can be cast in a suitable housing and functionality can be added before or after casting.

In accordance with the present invention, the structures of the present invention can be formed by a polymer phase inversion process, air casting (evaporation) and thermal inversion. For those systems with minimal or no adhesion, the formed structures can be separately prepared and inserted into the appropriate housing and held in place by mechanical means. In the preferred method, the formed structures are cast in situ in the desired housing. This results in the ability to include large amounts of media in the polymer matrix while still maintaining a three-dimensional porous structure. The membrane substructure serves as a support network enmeshing the particles, thus eliminating the need for frits or plugs, thereby

minimizing or even eliminating dead volume (the adsorptivity of the membrane may or may not participate in the adsorption process). However, porous frits plugs could be added if desired. Preferably the membranes or composite structures formed have an aspect ratio (average diameter to average thickness) of less than about 20, more preferably less than about 10, especially less than 1. For example, for adsorptive pipette tips, aspect ratios of two or less, more preferably less than 1 are preferred, especially between about 0.005-0.5. An aspect ratio within this range provides for suitable residence times of the sample in the composite structure during operation.

5

10

15

20

25

In the polymer phase inversion process, the solvent for the polymer must be miscible with the quench or inversion phase. For example, N-methyl-pyrolidone is a suitable solvent for polysulfones, polyethersulfones and polystyrene. In the latter case, polystryene pellets can be dissolved in N-methylpyrolidone and case-in-place. The resulting structure shows good adhesion to the walls of a polyolefin-based housing, and polysulfone. characteristics similar to adsorption Dimethylsulfoxide (DMSO), dimethylform-amide, butyrolactone, and N, N-dimethylacetamide sulfalane are also suitable solvents. (DMAC) is a suitable solvent for PVDF. Water is the preferred The polymer phase inversion process generally precipitant. results in an expansion of the structure to about two to three times its casting solution volume in the housing.

In the air casting process, a volatile solvent for the polymer binder is used. For example, in the case of cellulose acetate, acetone is a suitable volatile solvent. Air casting

generally results in a structure which is smaller than the casting solution volume. With this method, particles in the filled structures should be at least about  $30\mu$  to allow flow through the interstitial spaces after shrinkage without having to apply higher driving force.

5

10

15

20

25

The upper limit of particle amounts is dictated by casting solution viscosity. Depending on particle type, up to 40% (w/w) of particles can be added to the polymer without resulting in a casting solution too viscous to draw into the housing. Higher particle loadings may be achieved using higher temperature. Suitable particle sizes include particles in the range of from about 100 nanometers to about 100 microns in average diameter with or without porosity.

Suitable housing materials are not particularly limited, and include plastics (such as polyethylene and polypropylene), glass and stainless steel. Polyolefins, and particularly polypropylene, are preferred housing materials in view of the chemical adhesion that is created with the composite structure when the composite containing polysulfone, and in particular UDEL P3500 and P1700 polysulfones available from Amoco, is castin-place therein. Figure 16B illustrates such adhesion with a polypropylene pipette tip housing having a cast-in-place membrane therein prepared with spherical silica gel and polysulfone.

Suitable housing configurations are also not particularly limited, and include pipette tips, wells, multi-well arrays, plastic and glass cavities, sample preparation devices such as the MICROCON® microconcentrator, commercially available from

5

10

15

20

25

preferred housing The Corporation, etc.. Millipore configuration is substantially cylindrical, as the flow vectors during operation are substantially straight, similar chromatography, thereby minimizing or avoiding dilutional washing that might occur with non-cylindrical configurations. Although housings with volumes between about 0.1  $\mu$ l and about 5 mls. can be used for casting-in-place, volumes less than about 100 $\mu$ l are preferred, with volumes of from about 0.1-50 $\mu$ l, preferably from about  $0.2-20\mu l$ , are especially preferred. Pipette tip geometries having volumes as small as about 5 When chemical adhesion of the microliters can be used. composite structure to the housing walls is desired but is insignificant or non-existent, mechanical means can be used to maintain the composite structure in the housing. such as crimping, press fitting, heat shrinking the housing or a portion thereof, plasma treating the housing or a portion thereof, or chemically treating, such as etching, the housing or a portion thereof to promote adhesion. An advantage of adhesion to the housing wall is the ability to "seal" the composite structure Such sealing (by to the housing without mechanical means. whatever method) prevents the sample from channeling bypassing the composite during operation. the Preferably structures of the present invention have a final bed height of This allows for good washing, from about 0.05 to about 5mm. good density per unit volume, and results in a uniform precipitation during formation of the plug.

The structures of the present invention also can be castin-place in conventional multi-well arrays having suitable

Alternatively, as shown in Figures 5A-5D, multigeometries. well arrays 10 can be used as the housing, such as by casting the structures 11 of the present invention in place in the well Alternatively, Figure 5B shows an underdrain subassembly 13 having a plurality of wells 12 (enlarged in Figure 5D) with cast-in-place structures contained therein. The underdrain 13 can be adapted to be permanently or removably coupled to the reservoir array 10 by any suitable means, such as by snapping, so as to form removable "boot" assemblies containing the For convenience, each structures of the present invention. underdrain 13 can contain a polymer matrix having particles with different chemistry, so that the user chooses the appropriate underdrain 13 depending upon the application. Alternatively or in addition, the particle laden polymer matrix can differ from The reservoir housing 10 can be a plurality of well to well. open bores, or can include a membrane.

5

10

15

20

25

The composite structures and the micro sample preparation devices of the present invention containing the composite structures have a wide variety of applications, depending upon For example, applications include the particle selection. peptide and protein sample preparation prior to analysis, peptide removal from carbohydrate samples, amino acid clean-up immobilized enzymes for micro-volume analysis, prior to micro-affinity ligands for immobilized reactions, chromatography, isolation of supercoiled and cut plasmids, clean-up of PCR and DNA products, immobilized oligo dT for RNA isolation, dye terminator removal, sample preparation for Those skilled in the art will be able mental analysis. etc..

to choose the appropriate particles, polymer binder, particle chemistry and form geometry depending upon the desired application. In some cases, a mixture of particles can be used in the same devices. Alternatively or in addition, a multi-well device could have different chemistries for each separate well.

5

10

15

20

25

In the embodiment where the structures of the present invention are not filled with particles, symmetrical or asymmetrical semi-permeable structures, or a combination of symmetrical and asymmetrical semi-permeable structures, can be formed. In this embodiment, the preferred method of formation is casting in situ in the appropriate housing to form a self-retaining, self-supporting structure suitable for separations based on size or adsorption (depending on polymer identity). Functionality can be added to such a membrane to perform adsorption separations without the use of particles. For example, cellulose acetate can be treated with base to form cellulose, followed by an oxidant to render it reactive.

In the *in situ* formation process (either with filled or unfilled structures), the preferred method of formation involves precipitation by means of solvent exchange, such as by introducing the casting solution into the housing by any suitable means, such as where pressure is the driving force, for example by capillary action or by using a vacuum source. In the embodiment in which the housing is a pipette tip, a preferred driving force is a hand-held pipettor. Once the desired volume in the housing is filled with casting solution, the casting solution in the housing is contacted with a liquid in which the polymer is insoluble, preferably water, so that the polymer

precipitates in the housing. This can be accomplished by immersing the housing in the liquid, and/or drawing the liquid into the housing with a driving force such as by means of a vacuum. Through the exchange of water for the solvent, the structure precipitates. Those skilled in the art will appreciate that the solvent used to prepare the casting solution and the non-solvent can contain a variety of additives.

5

10

15

20

25

At the first contact of the polymer with the precipitant, there is virtually instaneous precipitation, thereby forming a semi-permeable barrier or "skin". Such a barrier is illustrated in Figure 17 as element 60 in a housing 62. This barrier slows the rate of further precipitation of the substructure. Once precipitation is complete, the initial semi-permeable barrier 60 can be removed, such as by cutting the housing at a point above the barrier at a point above the barrier or by abrading exposed polymer. The semi-permeable barrier 60 can be optionally left in place to carry out size-based separations with unfilled structures, as the barrier acts as a microfiltration membrane.

The cast in-place structure assumes the shape of the housing and results in a self-retaining homogeneous structure akin to a chromatographic column, providing a large surface area suitable for bind/elute chromatography (e.g., when particles are included in the polymer matrix) or for other analytical or biochemical techniques. Suitable driving forces include centrifugation, gravity, pressure or vacuum.

Without limitation, the following examples illustrate the objects and advantages of the present invention.

# Example 1: Strong Cation Exchange (SCX) Silica in 20 µl Pipette Tips

5

10

15

20

25 -

In a suitable small vessel, 5 grams of a 7 % (w/w) PVDF solution (Pennwalt Corp, KYNAR 761) was prepared in N,Ndimethyacetamide. To this, 1 gram of SCX, 200 Å, 15  $\mu$ m (Millipore, PN 85864) spherical silica was added and mixed thoroughly with a spatula. The mixture was allowed to equilibrate for 2 hours at room temperature, then mixed again. A 20  $\mu$ l fluted polypropylene disposable pipette tip was affixed to a common P-20 Pipetman (Gilson, Ranin, etc.) and the volume adjustment was set to 20  $\mu$ l. The plunger was depressed to the bottom and the end of the pipette was placed into the casting solution. While carefully watching, the plunger was slowly raised to fill the tip with ca. 0.5 - 1  $\mu$ l of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the pipette tip was removed and dipped into a bath of deionized water @ 60° C for ca. 5 seconds. After this brief period, pressure was released on the plunger and water was drawn into the tip to precipitate the polymer. When the water level was ca. 0.5 cm above the polymer height, the tip was ejected into the bath and solvent exchange was allowed to occur for ca. 5 minutes. The tip was removed from the water bath and any precipitated polymer located on the exterior was abraded off. The tip was re-affixed to the pipettor and the liquid expelled. flow is poor, ca. 0.25 mm can be cut off the end with a sharp To ensure that all solvent was removed, ca. 5 razor blade. to 20  $\mu$ l of deionized water was drawn in and expelled several times.

## Example 2: C18 Silica in Common 200 µl Pipette Tips

5

10

15

20

25

In a suitable small vessel, 5 grams of a 6% (w/w). polysulfone solution (Amoco, P3500) was prepared in N-methyl-2-pyrrolidone. To this 2 grams of C18, 200 Å, 15  $\mu m$ spherical silica (Millipore, PN 85058 ) was added and mixed thoroughly with a spatula. The mixture was allowed to equilibrate for 2 hours at RT., then mixed again. A 200  $\mu$ l fluted polypropylene disposable pipette tip was affixed to a common P-200 Pipetman (Gilson, Ranin, etc.) and the volume adjustment was set to 200  $\mu$ l. The plunger was depressed to the bottom and the end of the pipette was placed into the casting solution. While carefully watching, the plunger was slowly raised to fill the tip with ca. 2 - 5  $\mu$ l of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the tip was removed and dipped into a bath of deionized water at room temperature for ca. 5 seconds. After this brief period, pressure on the plunger was released and water was drawn into the tip to precipitate the polymer. When the water level was ca. 0.5 - 1 cm above the polymer height, the tip was ejected into the bath and solvent exchange was allowed to occur for ca. 5 minutes. The tip was removed from the water bath and any precipitated polymer located on the exterior was twisted off. The tip was re-affixed to the pipetter and the liquid expelled. flow is poor, ca. 0.5 mm can be cut off the end with a sharp razor blade. To ensure that all solvent was removed, ca. 50 to 200  $\mu$ l of deionized water was drawn in an expelled several

## Example 3: $60 \text{ Å}, 10 \mu \text{m}$ Normal Phase Silica in Wide Bore $1000\mu \text{l}$ Pipette Tips

5

10

15

20

25

In a suitable small vessel, 6 grams of a 6% (w/w) cellulose acetate solution (Eastman Kodak, 398-60) was prepared in N-methyl-2-pyrrolidone. To this, 1 gram of 60 Å, 10  $\mu\text{m}$  granular silica gel (Davison, Grade 710) was added and mixed thoroughly with a spatula. The mixture was allowed to equilibrate for 2 hours at room temperature, then mixed again. A wide bore 1000  $\mu$ l polypropylene pipette was affixed to a common P-1000 Pipetman (Gilson, Ranin, etc.) and the volume adjust was set to 1000  $\mu$ l. The plunger was depressed to the bottom and the end of the pipette was placed into the While carefully watching, the plunger was casting solution. slowly raised to fill the tip with ca. 10 - 25  $\mu l$  of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the tip was removed and dipped into a bath of deionized water for ca. 5 seconds. After this brief period, pressure on the plunger was released and water was drawn into the tip to precipitate the polymer. water level was ca. 1 cm above the polymer height, the tip was ejected into the bath and solvent exchange was allowed to take place for ca. 5 minutes. The tip was removed from the water bath and any precipitated polymer located on the exterior was abraded off. The tip was re-affixed to the pipettor and the liquid expelled. If the flow is poor, cut ca. 0.5 mm off the end with a sharp razor blade. To ensure that all solvent was removed, ca. 200 to 1000  $\mu l$  of deionized water was drawn in and expelled.

## Example 4: Fumed Silica in Wide Bore 200 $\mu$ l Pipette Tips

5

10

15

20

25

In a suitable small vessel, 8 grams of a 7.5 % (w/w) polysulfone solution (Amoco, P3500) was prepared in N-methyl-To this, 0.5 grams of fumed silica (Degussa, 2-pyrrolidone. Aerosil 200) were added and mixed thoroughly with a spatula. The mixture was allowed to equilibrate for 2 hours at room temperature, then mixed again. A 200  $\mu$ l wide bore polypropylene pipette was affixed to a common P-200 Pipetman (Gilson, Ranin, etc.) and the volume adjust was set to 200 The plunger was depressed to the bottom and the end of the pipette was placed into the casting solution. While carefully watching, the plunger was slowly raised to fill the tip with ca. 10 - 25  $\mu$ l of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the tip was removed and dipped into a bath of deionized water for ca. 5 seconds. After this brief period, pressure on the plunger was released and water was drawn into the tip to precipitate the polymer. When the water level was ca. 1 cm above the polymer height, the tip was ejected into the bath and solvent exchange was allowed to take place for ca. 5 The tip was removed from the water bath and any minutes. precipitated polymer located on the exterior was abraded off. The tip was re-affixed to the pipettor and the liquid expelled. If the flow is poor, cut ca. 0.5 mm off the end with a sharp razor blade. To ensure that all solvent was removed, ca. 200 to 1000  $\mu l$  of deionized water was drawn in and expelled.

# Example 5: C18, 15 um silica particle loaded membrane in place

5

10

15

20

25

In a small vessel, 5 grams of a 6 % (w/w) polysulfone solution (Amoco, P3500) was prepared in N-methyl-2pyrrolidone. To this, 2 grams of C18, 200 Å, 15  $\mu m$  silica (Millipore, PN 85864) was added and mixed thoroughly with a spatula. The mixture was allowed to equilibrate for 2 hours at room temperature, then mixed again. Using a pipette or eye dropper, 25 - 50  $\mu$ l of casting solution was dispensed into a suitable fixture. Examples of such devices include (but are not limited to) an Millipore Microcon or the wells of a 96 well filter plate. When preparing devices by this method, each chamber must contain a permeable barrier which will retain the solution (e.g. polypropylene fabric, membrane, etc.). Once added, the unit was gently tapped to ensure that the solution covered the entire barrier surface. The device was immersed in water for ca. 2 hours, and was gently stirred every 15 mins to promote solvent exchange. After this period, the units were removed and placed in either a centrifuge or vacuum manifold, as appropriate. cast in place structure was flushed with 500 to 1000  $\mu l$  of deionized water to ensure solvent removal.

# Example 6. Cast Porous End Plug in Wide Bore 1000 $\mu$ l Pipette Tips containing Loose 30 $\mu$ l Silica

In a suitable small vessel, 5 grams of a 7.5 % (w/w) polysulfone solution (Amoco, P3500) was prepared in N-methyl-2-pyrrolidone. The mixture was allowed to equilibrate for 2 hours at room temperature, then mixed again. A 1000  $\mu$ l wide

5

10

15

20

25

bore polypropylene pipette was affixed to a common P-1000 Pipetman (Gilson, Ranin, etc.) and the volume adjust was set to 1000  $\mu$ l. The plunger was depressed to the bottom and the end of the pipette was placed into the casting solution. While carefully watching, the plunger was slowly raised to fill the tip with ca. 2 - 10  $\mu l$  of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, the tip was removed and dipped into a bath of deionized water for ca. 5 seconds. After this brief period, pressure on the plunger was released and water drawn into the tip to precipitate the polymer. When the water level was ca. 0.5 cm above the polymer height, the tip was ejected into the bath and solvent exchange allowed to take place for ca. 5 The tip was removed from the water bath and any precipitated polymer located on the exterior was abraded off. The tip was re-affixed to the pipettor and the liquid expelled. If the flow is poor, cut ca. 0.5 mm off the end with a sharp razor blade. To ensure that all solvent was removed, ca. 100 to 500  $\mu$ l of deionized water was drawn in and expelled. The pipette was detached and any excess water in the upper chamber was removed with a cotton swab. mg of (250 Å) 30  $\mu m$  silica gel was weighed out and carefully added to the back end of the pipette. The pipette was tapped so that the silica rested on top of the cast-in-place If necessary, affix a suitable porous plug (cotton or polypropylene) in the upper chamber to prevent particle loss.

### Example 7: Cast Semi-Permeable Membrane Plug for Filtration

In a suitable vessel, 5 grams of 7.5% (w/w) polysulfone solution (Amoco, P3500) in N-methyl-2-pyrrolidone was The mixture is allowed to equilibrate for 2 hours prepared. at room temperature, and is then mixed again. A 1000  $\mu$ l wide bore polypropylene pipette is affixed to a common P-1000 Pipetman pipettor (Gilson, Ranin, etc.) and the volume adjust is set to 1000  $\mu$ l. The plunger is depressed to the bottom and the end of the pipette is placed into the casting solution. While carefully watching, the plunger was slowly raised to fill the tip with ca. 2-10  $\mu l$  of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the tip was removed, excess polymer solution was wiped off, and the tip was dipped into a bath of deionized water for about 5 seconds. After this brief period, pressure was released on the plunger and water was drawn into the tip to precipitate the polymer. When the water level was about 0.5 cm above the polymer height, the tip was ejected into the bath and solvent exchange was allowed to take place for about 5 minutes. The tip was reaffixed to the pipettor, the liquid expelled, and washed with 100-200  $\mu$ l of deionized water. When cast in this manner, the precipitated polymer had a semi-permeable skin at the orifice, which can be used as a filtration medium.

25

20

5

10

15

### Example 8: Porous 30 µl Silica Plugs Prepared by Evaporation

In a suitable vessel, 5 grams of a 10% (w/w) cellulose acetate solution (Eastman Kodak, 398-60) in acetone was

To this, 1 gram of methanol, 0.5 grams of prepared. deionized water and 1 gram of 250 Å, 30  $\mu m$  silica was added. The mixture was allowed to equilibrate for 2 hours at room temperature, and was then mixed again. A 1000  $\mu l$  wide bore polypropylene pipette was affixed to a common P-1000 Pipetman pipettor (Gilson) and the volume adjust was set to 1000  $\mu$ l. The plunger was depressed to the bottom and the end of the pipette was placed into the casting solution. The plunger was then slowly raised to fill the tip with about 5-10  $\mu$ l of casting solution. Once the tip contained sufficient liquid, equal pressure was maintained, and the tip was removed, excess fluid was wiped off, and the tip was placed in a rack to allow solvent to evaporate for about 16 hours. After this period, the tip was washed with about 10  $\mu l$  of distilled water.

5

10

15

# Example 9: 30 μl Silica End Plugs in Porous Polyethylene Prepared by Thermal Phase Inversion

In a suitable vessel, 5 grams of beaded polyethylene and 20 100 grams of mineral oil are added. The mixture is heated to 250°C on a hot plate with agitation. When the plastic liquifies, 4 grams of 250Å, 30  $\mu m$  silica is added and mixed thoroughly. Using a 1 ml graduated glass pipette with filler bulb, 50-100  $\mu$ l of the melt is drawn in. Once the tip 25 contains sufficient liquid, equal pressure is maintained, and the tip is removed, excess plastic is wiped off, the tip is The pipette is allowed to cool to room temperature. transferred to a methylene chloride bath for 1 hour to It is then removed, and the extract the mineral oil. 30

methylene chloride is expelled and allowed to air dry.

# Example 10. C18 Silica 200 µl Pipette Tips for Peptide Sample Preparation

Approximately 2.5  $\mu$ g of each peptide from a mixture 5 consisting of GlyTyr (1), ValTyrVal (2), methionine enkephalin (3), leucine enkaphalin (4) and angiotensin II (5) (in 100  $\mu$ l 0.1 % TFA) was adsorbed to a P200 pipette tip containing ca. 5  $\mu$ l of cast C18, 200Å, 15  $\mu$ m spherical silica. The solution was drawn up and expelled 4 times. 10 tip was then washed with 200  $\mu$ l of 0.1% TFA. Bound peptides were eluted with 80% Acetonitrile in 0.1% TFA/water. The eluted peptides were diluted with 4 parts of 0.1 % TFA and analyzed by reverse phase HPLC (linear acetonitrile gradient 5 - 30 % over 20 min). The resulting chromatogram was then 15 compared to that of the original mixture. (See Figures 6 and As expected, the GlyTyr, ValTyrVal, which are small and relatively hydrophilic, did not bind to the  $C_{18}$ . recoveries of the remaining 3 (adsorbed) peptides subsequent to elution ranged from 70-85 %. 20

# Example 11. Strong Cation Exchange Silica 200 $\mu$ l Pipette Tips

25

Approximately 2.5  $\mu$ g of each solute from a mixture consisting of a five peptides (see Example 10) (in 100  $\mu$ l in 10 % glacial acetic acid) were adsorbed to a P200 pipette tip containing ca. 5  $\mu$ l of cast, styrene sulfonate coated, 300 Å, 15  $\mu$ m spherical silica. Adsorption was performed during 4 complete uptake-withdraw cycles followed by a 100  $\mu$ l flush

with 20 % methanol/10 mM HCl. Bound sample was eluted with two 25  $\mu$ l volumes of 1.4 N ammonium hydroxide/50 % methanol. The eluted sample was analyzed by reversed phase HPLC and the resulting chromatogram was compared to that of the original mixture. (See Figures 6 and 8). The strong cation exchange tip bound all sample components, except GlyTyr. Such performance is consistent with the selectivity of sulfonic acid ion-exchange resins.

### Example 12. Immobilized Enzyme in 200 $\mu$ l Pipette Tips

5

10

15

20

25

Trypsin was covalently coupled to an aldehyde activated 300 Å, 15  $\mu$ m spherical silica and cast (20  $\mu$ l) into P200 tips for protein digestion in situ. Trypsin activity within the tip was assessed by monitoring the digestion of cytochrome via HPLC. A sample of cytochrome c (10  $\mu$ g in 100  $\mu$ l of 100 mM Tris, 1 mM CaCl<sub>2</sub> pH 8 @ 37C) was taken up into the tip for 15 minutes. The reaction was mixed 4X with a expel/draw cycle into an Eppendorf tube. The digest was analyzed by HPLC using a linear gradient of acetonitrile from 5 - 45% over 30 minutes (see Figure 10). The resulting chromatogram showed that greater than 90% of cytochrome c was digested after 15 minutes (see Figure 9 for undigested cytochrome c).

## Example 13: Immobilized Protein A in 200 µl Pipette Tips

Recombinant protein A was coupled to precast P200 tips containing aldehyde-activated 300 Å, 15  $\mu m$  spherical silica for the isolation of rabbit immunoglobulin (IgG). A 100  $\mu l$  sample of 1 mg/ml IgG and BSA in RIP buffer (150 mM NaCl, 1%

NP-40, 0.5% DOC, 0.1% SDS, 50 mM Tris, pH 8.0) was cycled six times through a tip containing 40  $\mu l$  of cast volume containing protein A immobilized beads. The tip was then washed with 5 volumes of RIP buffer prior to the elution. Desorption of bound IgG was performed with (two 25  $\mu$ l volumes) of 6M urea. The desorbed sample was diluted with 50  $\mu$ l of 2X SDS Laemmli sample buffer and boiled for 3 min prior to electrophoretic analysis. This protocol was also performed on a blank tip containing just polysulfone without beads which served as a background control. Electrophoresis was performed in a 10-16 % acrylamide gel shown (see Figure Samples are as follows: Lane 9: (MW marker); lanes 1-4: increasing amounts of protein A tip eluted sample; and lanes 5-8: increasing amounts of eluted IgG/BSA from the blank polysulfone tip. These results indicate selective binding of IgG to the Protein A tip with minimal nonspecific adsorption. Furthermore, the blank tip (lanes 5-8), in the presence of detergents (RIP buffer), did not exhibit adsorption of either IgG or BSA.

20

25

**15** 

5

10

## Example 14: 60 Å, 10 µm 1000 µl Pipette Tips for Supercoiled DNA

was grown in 3 - 5 ml of Luria broth containing plasmid pUC19 ampicillin at 37°C for 12-16 hours. 1.5 ml of the overnight culture was pelleted in a microfuge tube spun at maximum g-force for 30 sec at room temperature. Residual growth medium was removed while leaving the bacterial pellet intact. Plasmid DNA was then isolated using a modification of the

5

10

15

20

25

alkaline lysis procedure of Birnboim and Doly (Birnboim, H. C. and Doly, J. (1979). Nucleic Acids Res 7., 1513). Briefly, the bacterial pellet was resuspended by vortexing in 50  $\mu$ l of 50 mM glucose, 25 mM Tris-HCl (pH 8.0), 10 mM EDTA, and 10  $\mu$ g/ml RNase A. Next 100  $\mu$ l of 0.2 N NaOH, 1 % sodium dodecyl sulfate was added. The resulting suspension was incubated at room temperature for 2 min. Following the addition of 100  $\mu$ l of 3 M sodium acetate solution (pH 4.8), the suspension was mixed by vortexing then spun in a microfuge at maximum g-force for 2 min. The cleared lysate was transferred to a fresh microfuge tube to which 7 M guanidine hydrochloride (GuHCl) in 200 mM 2-(Nmorpholino) ethane sulfonic acid (MES) at pH 5.6 was added to a final concentration and volume of 4.4 M and 700  $\mu$ l, respectively. The resulting solution was drawn into a 1000  $\mu$ l polypropylene pipette tip with ca. 60  $\mu l$  of cast membrane containing ca. 60 Å, 10  $\mu m$  silica gel using a P-1000 pipettor. The solution was pipetted in-and-out for 2 - 2.5 minutes to allow extensive interaction between the DNA solution and the silica membrane matrix. The tip was then flushed once with 400  $\mu$ l of 80 % reagent grade alcohol. Residual alcohol is removed by repeated expulsion onto a paper towel. Plasmid DNA was eluted from the tip in 100  $\mu$ l of 10 mM Tris-HCl (pH 8.0), 1 mM EDTA (TE) by in-and-out pipetting 3x. Eluate fractions were adjusted to a final volume of 100  $\mu$ l with TE. Six tips were evaluated. quantitate plasmid DNA recovery, 20 % of the eluate, as well s 20% of the unbound filtrates, were analyzed by agarose gel

electrophoresis (See Figure 12). Included on the gel were samples of pUC19 plasmid DNA of known concentrations. (Lanes 1-4) Results of these experiments indicate that on average 2.5 mg of supercoiled plasmid was recovered (Lanes 5,7,9,11).

5

10

15

20

25

## Example 15: 60 Å, 10 µm Silica in Wide Bore 200 µl Pipette Tips for Linear DNA

The ability of 200  $\mu$ l polypropylene wide bore pipette tips containing ca. 20  $\mu$ l of cast 60 Å, 10  $\mu$ m silica-laden membrane to bind linearized DNA fragments (pBR322 digested with either BstNI or MspI, to generate DNA fragment ladders) or plasmid pBR322 DNA restricted with PstI and BamHI (generates large linear restriction fragments) was assessed. Five  $\mu g$  of linearized plasmid DNA was combined with GuHCl, pH 5.6 in MES to a final concentration of 0.5 M and volume of Prior to use, P-200 tips containing the silica membrane were pre-equilibrated in (2X) 200  $\mu$ l of 0.5 M GuHCl, The DNA/GuHCl solution was drawn into a pH 5.6 in MES. pipette tip and cycled in-and-out for 1.5-2.0 min to allow extensive interaction between the DNA binding mixture and the silica-laden membrane matrix. The tips were then washed with 125  $\mu$ l of 80% reagent grade alcohol to remove salts and other contaminants. Bound DNA was eluted from the tip matrix in 100  $\mu$ l TE, by in-and-out pipetting 3X. To measure DNA recovery, eluates and filtrates were analyzed by agarose gel electrophoresis (see Figure 13). In order to quantitate the amount of DNA recovered, samples representing 100%, 75%, 50%, and 25% of the starting material were run in lanes 1-4. Lanes 5, 7, 9, & 11 are the eluants. Estimate of band intensities

indicate recoveries in excess of 95%.

5

10

15

20

25

## Example 16. Fumed Silica in Wide Bore 200 µl Pipette Tips for PCR Amplified DNA

The ability of 200  $\mu$ l wide bore polypropylene pipette tips containing ca. 20  $\mu l$  of fumed silica immobilized in a polysulfone matrix was assessed for the purification of PCR amplified DNA (500 bp). Prior to use, tips were flushed 2X with 100  $\mu$ l of TE buffer and then equilibrated with 500  $\mu$ l of 3 M NaI in 200 mM MES buffer (pH 6.4). 50  $\mu$ l samples from the pooled PCR stock (ca. 3  $\mu g$  of DNA) were then combined with 7 M NaI to a final NaI concentration of 3.0 M. The total volume following addition of the NaI solution was 150  $\mu$ l. The sample was drawn in and expelled from the P-200 tips containing the cast fumed silica-laden membrane for 2-3 minutes allowing for extensive contact with the matrix. Each tip was then washed with 125  $\mu l$  of 80 % reagent grade alcohol to remove salts and other contaminants. Residual alcohol was removed by expelling the tip contents onto a paper towel. Bound PCR product was eluted in 50  $\mu$ l TE (pH 8.0). estimate DNA recovery, eluates and filtrates were analyzed by agarose gel electrophoresis (see Figure 14). Loads representing 100%, 75%, 50%, and 25% of the starting material were run in lanes 1-4 as controls. Note the presence of the lower band which indicates a slight primer-dimer contamination. The use of immobilized fumed silica along with NaI appears to give an amplified DNA recovery in excess of 90 %. In addition, there appears to be a reduction in the primer-dimer contaminant. (See Lanes 5,7,9,11).

# Example 17. Cast Porous End Plug with Loose 30 Micron Silica in a 200 µl Pipette Tip for DNA Isolation

5

10

15

20

200  $\mu$ l pipette tips containing ca 5 - 10  $\mu$ l of cast (7.5 %) polysulfone as a porous end plug and 2 - 4 mg of loose 250 Å, 30  $\mu m$  silica was assayed for the ability to bind linear and supercoiled plasmid DNA. Regarding linear DNA, approximately 5  $\mu g$  of pBR322 was first digested with MspI in 45  $\mu$ l TE (10 mM Tris-HCl, 1 mM EDTA), pH 8.0, and then combined with 100  $\mu$ l of 7 M guanidine hydrochloride (GuHCl) in 200 mM MES buffer at pH 5.6. The final concentration of GuHCl in the solution was 4.7 M. The resulting solution was drawn (once) into a 200  $\mu$ l pipette tip and allowed to extensively contact the silica by inverting the pipetman with the affixed tip for approximately 2 min. The DNA adsorbed to the tips was then washed and eluted as described in Example Loads representing 100%, 75%, 50% and 25% of the starting material where run in Lanes 1-4 as controls. Results from experiments using this format indicate that DNA recoveries of better than 75% can be achieved (see Figure 15, Lanes 5 and 7).

What is claimed is:

5

10

15

20

25

1. A housing defining a volume, said housing containing in a portion of said volume a three dimensional structure comprising a plurality of sorptive particles entrapped in a porous polymer matrix, said structure having an aspect ratio of less than about 10.

- 2. The housing of claim 1, wherein said housing has a first open end and a second open end spaced from said first open end, and wherein said three dimensional structure is contiguous with said second open end.
- 3. The housing of claim 1, wherein said aspect ratio is less than 5.
- 4. The housing of claim 1, wherein said portion of said volume containing said structure is from about 0.1 microliters to about 10 milliliters.
- 5. The housing of claim 1, wherein said housing is a pipette tip having a first end and a second end spaced from said first end with said structure in between.
- 6. The housing of claim 5, wherein said first end has an internal diameter larger than the internal diameter of said second end.
  - 7. The housing of claim 1, wherein said polymer is selected from the group consisting of polysulfone, polyethersulfone, polytetrafluoroethylene, cellulose acetate, polystyrene, polystyrene/acrylonitrile copolymer and PVDF.
  - 8. The housing of claim 7, wherein said polymer is polysulfone, and wherein said housing comprises a polyolefin.

9. The housing of claim 8, wherein said polyolefin is polypropylene.

- 10. The housing of claim 1, wherein said particles are selected from the group consisting of polystryrenedivinylbenzene beads, functionalized polystryrenedivinylbenzene beads, silica, fumed silica, derivitized silica and activated carbon.
- 11. The housing of claim 1, wherein said three dimensional structure comprises a semi-permeable barrier.
- 12. A method of casting a composite structure in a housing, said method comprising:

forming a solution of a polymer;

5

15

20

adding particles to said solution to form a casting solution;

introducing said casting solution into said housing; and subjecting said casting solution to a phase inversion so as to cause said polymer to form a porous polymer matrix comprising said particles.

- 13. The method of claim 12, wherein said solution is formed by dissolving said polymer in a solvent.
- 14. The method of claim 12, wherein said solution is formed by dissolving said polymer in a mixture of a solvent and non-solvent for said polymer.
- 15. The method of claim 12, wherein said phase inversion is caused by contacting said casting solution in said housing with a liquid in which said polymer is insoluble.

16. The method of claim 12, wherein said phase inversion is caused by evaporation of said solvent.

- 17. The method of claim 12, further comprising removing said porous polymer matrix from said housing and introducing said porous polymer matrix into a second housing.
- 18. A method of casting a membrane in a housing, said method comprising:

forming a solution of a polymer;

5

10

15

20

25

introducing said solution into said housing; and subjecting said solution to a phase inversion so as to cause said polymer to precipitate in said housing.

- 19. The method of claim 18 wherein said solution is formed by dissolving said polymer in a solvent.
- 20. The method of claim 18, wherein said phase inversion is caused by contacting said solution in said housing with a liquid in which said polymer is insoluble.
- 21. The method of claim 18, further comprising removing any semi-permeable barrier that is formed during said precipitation.
- 22. A housing defining a volume, said housing containing in a portion of said volume a three dimensional structure comprising a porous polymer matrix, said structure having an aspect ratio of less than about 10.
- 23. The housing of claim 22, wherein said housing has a first open end and a second open end spaced from said first open end, and wherein said structure is contiguous with said second open end.

WO 98/37949 PCT/US98/03696

24. The housing of claim 22, wherein said aspect ratio is less than 5.

25. The housing of claim 22, wherein said portion of said volume containing said structure is from about 0.1 microliters to about 10 milliliters.


5


10

15

- 26. The housing of claim 22, wherein said housing is a pipette tip having a first end and a second end spaced from said first end with said structure in between.
- 27. The housing of claim 26, wherein said first end has an internal diameter larger than the internal diameter of said second end.
- 28. The housing of claim 22, wherein said polymer is selected from the group consisting of polysulfone, polyethersulfone, polytetrafluoroethylene, cellulose acetate, polystyrene, polystyrene/acrylonitrile copolymer and PVDF.
- 29. The housing of claim 28, wherein said polymer is polysulfone, and wherein said housing comprises a polyolefin.
- 30. The housing of claim 29, wherein said polyolefin is polypropylene.
- 20 31. The housing of claim 22, wherein said three dimensional structure comprises a semi-permeable barrier.

1/18








FIG. 2B




FIG. 3A




FIG. 3B




FIG. 4





FIG. 5A






FIG. 5C






FIG. 5D












FIG. 8







# 1 2 3 4 5 6 7 8 9 10 11 12



## p. 98 (3885)


#### 1ane

- 1. 100% recovery equivalent
- 2. 75% recovery equivalent
- 3. 50% recovery equivalent
- 4. 25% recovery equivalent
- 5. Tip #1 eluate
- 6. Tip #1 flow through

## lane

- 7. Tip #2 eluate
- 8. Tip #2 flow through
- 9. Tip #3 eluate
- 10. Tip #3 flow through
- 11. Tip #4 eluate
- 12. Tip #4 flow through

FIG. 12



# 1 2 3 4 5 6 7 8 9 10 11 12


## lane

- 1. 100% recovery equivalent
- 2. 75% recovery equivalent
- 3. 50% recovery equivalent
- 4. 25% recovery equivalent
- 5. Tip #1 eluate
- 6. Tip #1 flow through

#### lane

- 7. Tip #2 eluate
- 8. Tip #2 flow through
- 9. Tip #3 eluate
- 10. Tip #3 flow through
- 11. Tip #4 eluate
- 12. Tip #4 flow through

FIG. 13



1 2 3 4 5 6 7 8 9 10 11 12

## 1ane

- 1. 100% recovery equivalent
- 2. 75% recovery equivalent
- 3. 50% recovery equivalent
- 4. 25% recovery equivalent
- 5. Tip #1 eluate
- 6. Tip #1 flow through

#### 1ane

- 7. Tip #2 eluate
- 8. Tip #2 flow through
- 9. Tip #3 eluate
- 10. Tip #3 flow through
- 11. Tip #4 eluate
- 12. Tip #4 flow through

FIG. 14

1 2 3 4 5 6 7 8




FIG. 15

#### LANE

- 1. 100% recovery equivalent
- 2. 75% recovery equivalent
- 3. 50% recovery equivalent
- 4. 25% recovery equivalent
- 5. Tip #1 eluate
- 6. Tip #1 flow through
- 7. Tip #2 eluate
- 8. Tip #2 flow through

93X

700µM




FIG. 16A

91X

720µM




FIG. 16B

-62

380X




FIG. 17

# INTERNATIONAL SEARCH REPORT

International application No. PCT/US98/03696

A. CLASSIFICATION OF SUBJECT MATTER			
IPC(6) :B01D 61/00, 63/00			
US CL :210/321.6, 321.72, 321.75, 502.1, 500.41, 500.36, 500.35, 500.42,			
According to International Patent Classification (IPC) or to both national classification and IPC			
B. FIELDS SEARCHED			
Minimum documentation searched (classification system followed by classification symbols)			
U.S. : 210/321.6, 321.72, 321.75, 502.1, 500.41, 500.36, 500.35, 500.42; 264/41			
Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched			
Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)			
C. DOCUMENTS CONSIDERED TO BE RELEVANT			
Category*	Citation of document, with indication, where ap	propriate, of the relevant passages	Relevant to claim No.
Y	US 5,556,598 A (RAYBUCK ET AL 96), entire disclosure.	) 17 September 1996 (17-09-	1-31
Y	US 5,582,892 A (ANDERSON) 10 December 1996 (10-12-96), 1-31 entire disclsorure.		
Y	US 4,810,381 A (HAGEN ET AL) 07 March 1989 (07-03-89), 7-11, 28-31 entire disclosure.		
Y	US 5,464,541 A (AYSTA ET AL) 07 November 1995 (07-11-95), entire disclosure.		
Y,P	US 5,637,506 A (GOKEN ET AL) 10 disclosure.	1-31	
Y	US 5,558,771 A (HAGEN ET AL) 24 September 1996 (24-09-96), 1-31 entire disclosure.		
Further documents are listed in the continuation of Box C. See patent family annex.			
* Special categories of cited documents:  "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention			
to	be of particular relevance	"X" document of particular relevance; the	
*L* do	tier document published on or after the international filing date nument which may throw doubts on priority claim(s) or which is ad to establish the publication date of another citation or other	considered novel or cannot be consider when the document is taken alone	
special reason (as specified)  "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is			
*O" document referring to an oral disclosure, use, exhibition or other combined with one or more other such documents, such combination being obvious to a person skilled in the art			
	cument published prior to the international filing data but later than priority data claimed	"A" document member of the same patent	family
Date of the actual completion of the international search  Date of mailing of the international search report			
25 MARCH 1998 0 4 JUN 1998			
Name and mailing address of the ISA/US Commissioner of Patents and Trademarks Box PCT Washington, D.C. 20231  Authorized officer ANA FORTUNA			
Facsimile No. (703) 305-3230 Telephone No. (703) 308-2351			